




# VERONIKIS INTRA-FIBEROPTIC PELVIC ASSISTANT

**The Most Exciting New  
Instrument for Use in:**

- Urethropexy
- Sacrocolpopexy
- Paravaginal Repair
- Vaginal Surgery
- Rectovaginal Septum  
Endometriosis




VIPA comes complete with fiber optic cable, 4 VIPA tips.

In Urethropexy Procedures and Paravaginal Repairs: The VIPA provides the guiding light for accurate suture placement by illuminating the anatomy lateral to the bladder neck.



VIPA Urethropexy Tip

In Sacrocolpopexy: The VIPA steadies the vaginal vault and provides a suture platform with fiberoptic illumination.




VIPA Vaginal Cuff Tips

In Vaginal Surgery: The VIPA illuminates the rectum to visualize endometriosis lesions on the rectovaginal septum.



VIPA  
Rectal  
Probe



**No Charge  
30-Day Evaluation  
Available!**

## **SUPERIOR MEDICAL LIMITED**

*Canada's OB/GYN Company*

*520 Champagne Drive*

*Toronto, Ontario*

*Canada M3J 2T9*

Telephone: (416) 635-9797

Toll Free: 1-800-268-7944

Fax: (416) 635-8931

e-mail: [info@superiormedical.com](mailto:info@superiormedical.com)

Web site: [www.superiormedical.com](http://www.superiormedical.com)

Product Description	Catalog No.
VIPA	.350-370
VIPA Urethropexy Replacement Tip	.350-374
VIPA Vaginal Cuff Tips (small)	.350-376
VIPA Vaginal Cuff Tips (large)	.350-378
VIPA Rectal Probe	.350-380
VIPA Fiber Optic Cable - (Storz to ACMI Connection)	.350-372